REFORPOST®

PORTEGUÊS

INDICAÇÕES

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ENGLISH

INDICATIONS

For the support of restorations and prosthetic crowns.

COMPOSITION

Peeso drills.....	Stainless steel
Rubber Stops.....	Latex
Template for radiographic measurement.....	Plastic
REFORPOST® GLASS FIBER:	Glass fiber (80%); pigmented resin (19%); stainless steel filament (1%).
REFORPOST® CARBON FIBER:	Carbon fiber (72%); pigmented resin (22%); stainless steel filament (6%).

TECHNIQUE OF USE

A. ROOT CANAL PREPARATION

1. Select the pin using the template and the periapical x-ray of the tooth (photo 1);
2. Remove some of the root canal filling with heated instruments to the predetermined depth. The use of drills can cause apical displacement of the gutta-percha. Allow at least 3 to 5 mm of filling at the root apex due to the presence of accessory canals, which must remain sealed to prevent bacterial contamination (photo 2);
3. Prepare the canal with the drill corresponding to the diameter of the selected post (photo 3).

B. PREPARATION OF REFORPOST FOR CEMENTATION

1. Position the pin to verify its adaptation. Make a mark on the pin 2 to 3 mm from the antagonist tooth (photo 4) - and cut it with a diamond bit at high rotation under cooling. Do not use scissors or pliers, because the maceration of the fibers changes the mechanical properties of the pin;
2. Clean the post with alcohol and dry it (photo 5);
3. Apply a layer of silane (SILANO ANGELUS®) and let it dry for 1 minute (photo 6);
- ④ WARNING: It is not necessary to apply the silane on REFORPOST® CARBON FIBER.
4. Apply the ANGELUS® FUSION-DURALINK® CATALYST to the post according to its instructions (photo 7).

C. CEMENTATION

1. Clean the canal with the air-water spray and dry it with air and absorbent paper points (photo 8);
2. Etch the canal with ANGELUS® 37% PHOSPHORIC ACID for 15 seconds (photo 9);
3. Wash again the canal and dry it with air and absorbent paper points (photo 10). Apply ANGELUS® FUSION-DURALINK® PRIMER in the canal and wait for 1 minute;
4. Apply ANGELUS® FUSION-DURALINK® CATALYST according to the instructions. The adhesive components must not be mixed before application. They are polymerized chemically without light curing (photo 11);
5. Spatulate the resinous cement selected* and fill the conduit;
6. Put the post in the canal, remove excess of the cement and wait for the polymerization. If the cement is dual, light cure it (photo 13).

D. CROWN BUILD-UP

1. Do another etching with ANGELUS® 37% PHOSPHORIC ACID and apply an adhesive on the remaining tooth structures and post. Build the crown core with composite resin (photos 14 and 15), or with ANGELUS® REFORCORE®. REFORCORE® is a prefabricated glass fiber core, for the technique of direct cores (photos 16 and 17).

E. USE OF ACCESSORY POSTS IN WEAKENED ROOTS AND WIDE CANALS

1. Fill the remaining spaces of the canal with ANGELUS® REFORPIN® accessory posts. The decrease in thickness of the cement and the better retention provided by REFORPIN® reduce the risk of post/core loosening (photos 18 and 19).

WARNINGS | PRECAUTIONS

- Use of direct cores with fiber posts is not indicated for extensive oral rehabilitation in patients with parafunctional habits;
 - Use of the direct core with fiber post should only be indicated to teeth with a minimum height of 2 mm of coronal dentin;
 - Fiber posts should only be cut with water refrigerated high speed diamond burs. Do not use scissors or cutting pliers. Use goggles and a mask to prevent inhalation of particles generated during cutting;
 - Fiber posts must be completely covered by the cementing agent or the composite resin. Exposed fibers absorb water which deteriorates the mechanical properties of the posts.
- ① IMPORTANT: The indication of any post depends on careful evaluation of the quality of endodontic treatment and periapical health.

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